



Product description

Order number	72.688.004
Product description	Reaction tube, working volume: 1.5 ml, material: PP, transparent, cap: yellow, retaining loop cap, cap attached, 500 piece(s)/bag

Product characteristics

Graduation	no
Cap	retaining loop cap
Volume of work	1,5 ml
Centrifugation max (RCF)	20000 x g

Size

Diameter	10.8 mm
Length of product	39 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	yellow

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	5000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	5000
Piece(s) / pallet	180000
Depth of case	398 mm
Width of case	398 mm
Height of case	300 mm
Case volume	0.0475 cbm
Weight of product	0.0012 kg
Weight of case	6.59 kg
EAN of inner box	4038917116404
EAN case	4038917015097

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

