



### **Product description**

|                     |   |
|---------------------|---|
| Order number        | 72.689.001  |
| Product description | Reaction tube, working volume: 2 ml, material: PP, transparent, cap: blue, retaining loop cap, cap attached, 500 piece(s)/bag |

### **Product characteristics**

|                          |                    |
|--------------------------|--------------------|
| Graduation               | no                 |
| Cap                      | retaining loop cap |
| Volume of work           | 2 ml               |
| Centrifugation max (RCF) | 20000 x g          |

### **Size**

|                   |         |
|-------------------|---------|
| Diameter          | 10.8 mm |
| Length of product | 39 mm   |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

|                   |                                  |
|-------------------|----------------------------------|
| Product material  | Polypropylene (PP)               |
| Colour of product | transparent                      |
| Closure material  | Low Density Polyethylene (LD-PE) |
| Colour of cap     | blue                             |

### **Purity & certification**

|                  |                               |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified     | CE                            |
| Batched          | yes                           |

### **Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 5000          |
| Type of smallest subpackaging | bag           |
| Piece(s) / inner box          | 500           |
| Piece (s) / outer case        | 5000          |
| Piece(s) / pallet             | 150000        |
| Depth of case                 | 398 mm        |
| Width of case                 | 398 mm        |
| Height of case                | 340 mm        |
| Case volume                   | 0.0539 cbm    |
| Weight of product             | 0.0014 kg     |
| Weight of case                | 7.48 kg       |
| EAN of inner box              | 4038917116428 |
| EAN case                      | 4038917015110 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

