



Product description

Order number	72.690.301
Product description	Reaction tube, working volume: 1.5 ml, material: PP, transparent, cap: transparent, cap attached, with injected graduation and writing space, 500 piece(s)/bag

Product characteristics

Graduation	yes
Volume of work	1,5 ml
Centrifugation max (RCF)	20000 x g
Autoclavable (max. °C)	121 °C

Size

Sample volume	1.5 ml
Length of product	39 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Reaction tube, 1.5 ml, PP

Page 2

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Colour of cap	transparent

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	5000
Piece(s) / pallet	160000
Case volume	0.0558 cbm
Weight of product	0.0011 kg
Weight of case	5.96 kg
EAN of inner box	4038917369268
EAN case	4038917369251

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-26

This document was prepared by EDP support and is valid without signature.



SARSTEDT