



Product description

Order number	72.690.304
Product description	Reaction tube, working volume: 1.5 ml, material: PP, brown, cap: brown, cap attached, with injected graduation and writing space, 500 piece(s)/bag

Product characteristics

Graduation	yes
Volume of work	1,5 ml
Centrifugation max (RCF)	20000 x g
Autoclavable (max. °C)	121 °C

Size

Sample volume	1.5 ml
Length of product	39 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Reaction tube, 1.5 ml, PP

Page 2

Material & colours

Product material	Polypropylene (PP)
Colour of product	brown
Colour of cap	brown

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	5000
Piece(s) / pallet	160000
Case volume	0.0558 cbm
Weight of product	0.001 kg
Weight of case	5.6 kg
EAN of inner box	4038917369299
EAN case	4038917369282

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2024-03-28

This document was prepared by EDP support and is valid without signature.



SARSTEDT