



### **Product description**

Order number	72.690.375
Product description	Reaction tube, working volume: 1.5 ml, material: PP, blue, cap: blue, cap attached, with injected graduation and writing space, 500 piece(s)/bag

### **Product characteristics**

Graduation	yes
Volume of work	1,5 ml
Centrifugation max (RCF)	20000 x g
Autoclavable (max. °C)	121 °C

### **Size**

Sample volume	1.5 ml
Length of product	39 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Product Specification**  
**Reaction tube, 1.5 ml, PP**

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	blue
Colour of cap	blue

**Purity & certification**

Product category	no medical device   IVD
Batched	yes

**Packaging**

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	5000
Piece(s) / pallet	160000
Case volume	0.0558 cbm
Weight of product	0.001 kg
Weight of case	5.6 kg
EAN of inner box	4038917369442
EAN case	4038917369435

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

