



Product description

| | |
|---------------------|--|
| Order number | 72.690.999 |
| Product description | Reaction tube, working volume: 1.5 ml, material: PP, colour mix, cap: colour mix, cap attached, with injected graduation and writing space, 500 piece(s)/bag |

Product characteristics

| | |
|--------------------------|-----------|
| Graduation | yes |
| Volume of work | 1,5 ml |
| Centrifugation max (RCF) | 20000 x g |
| Autoclavable (max. °C) | 121 °C |

Size

| | |
|-------------------|--------|
| Sample volume | 1.5 ml |
| Length of product | 39 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

| | |
|-------------------|--------------------|
| Product material | Polypropylene (PP) |
| Colour of product | colour mix |
| Colour of cap | colour mix |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1000 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 100 |
| Piece (s) / outer case | 5000 |
| Piece(s) / pallet | 140000 |
| Depth of case | 410 mm |
| Width of case | 461 mm |
| Height of case | 345 mm |
| Case volume | 0.0649 cbm |
| Weight of case | 5.19 kg |
| EAN of inner box | 4038917116756 |
| EAN case | 4038917015288 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

