



Product description

| | |
|---------------------|---|
| Order number | 72.695.002 |
| Product description | SafeSeal reaction tube, working volume: 2 ml, material: PP, red, cap: red, SafeSeal cap, cap attached, with injected graduation and writing space, 250 piece(s)/bag |

Product characteristics

| | |
|--------------------------|--------------|
| Graduation | yes |
| Cap | SafeSeal cap |
| Volume of work | 2 ml |
| Centrifugation max (RCF) | 25000 x g |

Size

| | |
|-------------------|---------|
| Diameter | 10.8 mm |
| Length of product | 40 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

| | |
|-------------------|--------------------|
| Product material | Polypropylene (PP) |
| Colour of product | red |
| Closure material | Polypropylene (PP) |
| Colour of cap | red |

Purity & certification

| | |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified | CE |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 2000 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 1000 |
| Piece (s) / outer case | 2000 |
| Piece(s) / pallet | 100000 |
| Depth of box | 301 mm |
| Width of box | 221 mm |
| Height of box | 161 mm |
| Depth of case | 458 mm |
| Width of case | 313 mm |
| Height of case | 180 mm |
| Case volume | 0.0258 cbm |
| Weight of product | 0.0011 kg |
| Weight of case | 3.09 kg |
| EAN of inner box | 4038917262750 |
| EAN case | 4038917262682 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

