



Product Specification
SafeSeal reaction tube, 2 ml, PP, PCR
Performance Tested

Page 1



Product description

Order number	72.695.400
Product description	SafeSeal reaction tube, working volume: 2 ml, material: PP, transparent, cap: transparent, SafeSeal cap, cap attached, with injected graduation and writing space, PCR Performance Tested, 250 piece(s)/minigrip bag

Product characteristics

Graduation	yes
Cap	SafeSeal cap
Volume of work	2 ml
Centrifugation max (RCF)	25000 x g

Size

Diameter	10.8 mm
Length of product	40 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
SafeSeal reaction tube, 2 ml, PP, PCR
Performance Tested

Page 2

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	Polypropylene (PP)
Colour of cap	transparent

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Minimum order qty.	2000
Type of smallest subpackaging	minigrip bag
Piece(s) / inner box	1000
Piece (s) / outer case	2000
Piece(s) / pallet	100000
Depth of box	301 mm
Width of box	221 mm
Height of box	161 mm
Depth of case	458 mm
Width of case	313 mm
Height of case	180 mm
Case volume	0.0258 cbm
Weight of product	0.0011 kg
Weight of case	3.15 kg
EAN of inner box	4038917223188
EAN case	4038917223171

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-20

This document was prepared by EDP support and is valid without signature.



SARSTEDT