



### **Product description**

Order number	72.696
Product description	Reaction tube, working volume: 1.5 ml, material: PP, transparent, without cap, 500 piece(s)/bag

### **Product characteristics**

Graduation	no
Volume of work	1,5 ml
Centrifugation max (RCF)	20000 x g

### **Size**

Diameter	10.8 mm
Length of product	39 mm

### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Reaction tube, 1.5 ml, PP**

**Page 2**

**Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

**Packaging**

Minimum order qty.	10000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	10000
Piece(s) / pallet	240000
Depth of case	398 mm
Width of case	398 mm
Height of case	415 mm
Case volume	0.0657 cbm
Weight of product	0.0009 kg
Weight of case	9.43 kg
EAN of inner box	4038917083485
EAN case	4038917002400

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-26

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**