



Product Specification
SafeSeal reaction tube, 5 ml, PP, PCR
Performance Tested

Page 1

Product description

Order number	72.701.400
Product description	SafeSeal reaction tube, working volume: 5 ml, material: PP, transparent, cap: natural, SafeSeal cap, cap attached, with injected graduation and writing space, PCR Performance Tested, 100 piece(s)/minigrip bag

Product characteristics

Graduation	yes
Cap	SafeSeal cap
Volume of work	5 ml
Centrifugation max (RCF)	25000 x g

Size

Diameter	17 mm
Length of product	60 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	Polypropylene (PP)
Colour of cap	natural

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Minimum order qty.	200
Type of smallest subpackaging	minigrip bag
Piece(s) / inner box	200
Piece (s) / outer case	800
Piece(s) / pallet	38400
Depth of box	294 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
SafeSeal reaction tube, 5 ml, PP, PCR
Performance Tested

Width of box	184 mm
Height of box	124 mm
Depth of case	598 mm
Width of case	398 mm
Height of case	145 mm
Case volume	0.0345 cbm
Weight of product	0.0031 kg
Weight of case	3.46 kg
EAN of inner box	4038917398558
EAN case	4038917398541

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

