



### **Product description**

Order number	72.609.311
Product description	Screw cap micro tube, working volume: 2 ml, skirted conical base, with knurling, transparent, without cap, sterile, 500 piece(s)/double bag

### **Product characteristics**

Graduation	no
Knurls	yes
Base shape	skirted conical base
Volume of work	2 ml
Centrifugation max (RCF)	20000 x g

### **Size**

Diameter	10.8 mm
Length of product	44 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Screw cap micro tube, 2 ml, sterile**

**Page 2**

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

**Purity & certification**

Satisfies the requirement	IATA
Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

**Packaging**

Minimum order qty.	500
Type of smallest subpackaging	double bag
Piece(s) / inner box	500
Piece (s) / outer case	2500
Piece(s) / pallet	140000
Case volume	0.029 cbm
Weight of product	0.0012 kg
Weight of case	3.63 kg
EAN of inner box	4038917434379
EAN case	4038917434362

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-05

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**