



Product description

Order number	72.687.772
Product description	Screw cap micro tube, working volume: 1.5 ml, conical base, transparent, without cap, sterile, 500 piece(s)/double bag

Product characteristics

Graduation	no
Knurls	no
Base shape	conical base
Volume of work	1,5 ml
Centrifugation max (RCF)	20000 x g

Size

Diameter	10.8 mm
Length of product	44 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Satisfies the requirement	IATA
Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order qty.	2500
Type of smallest subpackaging	double bag
Piece(s) / inner box	500
Piece (s) / outer case	2500
Piece(s) / pallet	130000
Depth of case	548 mm
Width of case	368 mm
Height of case	130 mm
Case volume	0.0262 cbm
Weight of product	0.0011 kg
Weight of case	3.27 kg
EAN of inner box	4038917116329
EAN case	4038917043823

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

