



### **Product description**

|                     |                                                                                                                                                         |
|---------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------|
| Order number        | 72.692.105                                                                                                                                              |
| Product description | Screw cap micro tube, working volume: 1.5 ml, conical base, with knurling, transparent, cap: natural, cap attached assembled, sterile, 100 piece(s)/bag |

### **Product characteristics**

|                          |                    |
|--------------------------|--------------------|
| Graduation               | no                 |
| Cap                      | retaining loop cap |
| Closure type             | screw cap          |
| Knurls                   | yes                |
| Base shape               | conical base       |
| Volume of work           | 1,5 ml             |
| Centrifugation max (RCF) | 20000 x g          |

### **Size**

|                   |         |
|-------------------|---------|
| Diameter          | 10.8 mm |
| Length of product | 44 mm   |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

|                   |                                        |
|-------------------|----------------------------------------|
| Product material  | Polypropylene (PP)                     |
| Colour of product | transparent                            |
| Closure material  | Polypropylene (PP)                     |
| Colour of cap     | natural                                |
| O-ring material   | Ethylene propylene diene rubber (EPDM) |

### **Purity & certification**

|                           |                               |
|---------------------------|-------------------------------|
| Satisfies the requirement | IATA                          |
| Product category          | In vitrodiagnostic device, CE |
| CE certified              | CE                            |
| Purity standard           | sterile                       |
| Sterilisation             | Electron irradiation          |
| Batched                   | yes                           |

### **Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 1000          |
| Type of smallest subpackaging | bag           |
| Piece(s) / inner box          | 100           |
| Piece (s) / outer case        | 1000          |
| Piece(s) / pallet             | 96000         |
| Depth of case                 | 348 mm        |
| Width of case                 | 278 mm        |
| Height of case                | 140 mm        |
| Case volume                   | 0.0135 cbm    |
| Weight of product             | 0.0015 kg     |
| Weight of case                | 1.81 kg       |
| EAN of inner box              | 4038917116879 |
| EAN case                      | 4038917015387 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

