



Product description

| | |
|---------------------|--|
| Order number | 72.730.100 |
| Product description | Screw cap micro tube, working volume: 0.5 ml, skirted conical base, with knurling, transparent, cap: natural, cap attached, 500 piece(s)/bag |

Product characteristics

| | |
|--------------------------|----------------------|
| Graduation | no |
| Cap | retaining loop cap |
| Knurls | yes |
| Base shape | skirted conical base |
| Volume of work | 0,5 ml |
| Centrifugation max (RCF) | 20000 x g |

Size

| | |
|-------------------|---------|
| Diameter | 10.8 mm |
| Length of product | 44 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

| | |
|-------------------|--|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |
| Closure material | Polypropylene (PP) |
| Colour of cap | natural |
| O-ring material | Ethylene propylene diene rubber (EPDM) |

Purity & certification

| | |
|---------------------------|-------------------------------|
| Satisfies the requirement | IATA |
| Product category | In vitrodiagnostic device, CE |
| CE certified | CE |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 5000 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 500 |
| Piece (s) / outer case | 5000 |
| Piece(s) / pallet | 120000 |
| Depth of case | 398 mm |
| Width of case | 398 mm |
| Height of case | 415 mm |
| Case volume | 0.0657 cbm |
| Weight of product | 0.0016 kg |
| Weight of case | 8.7 kg |
| EAN of inner box | 4038917119450 |
| EAN case | 4038917016643 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

