



### Product description

Order number	72.730.105
Product description	Screw cap micro tube, working volume: 0.5 ml, skirted conical base, with knurling, transparent, cap: natural, cap attached assembled, sterile, 100 piece(s)/bag

### Product characteristics

Graduation	no
Cap	retaining loop cap
Closure type	screw cap
Knurls	yes
Base shape	skirted conical base
Volume of work	0,5 ml
Centrifugation max (RCF)	20000 x g

### Size

Diameter	10.8 mm
Length of product	44 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	Polypropylene (PP)
Colour of cap	natural
O-ring material	Ethylene propylene diene rubber (EPDM)

### **Purity & certification**

Satisfies the requirement	IATA
Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

### **Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	96000
Depth of case	348 mm
Width of case	278 mm
Height of case	140 mm
Case volume	0.0135 cbm
Weight of product	0.0016 kg
Weight of case	1.9 kg
EAN of inner box	4038917119474
EAN case	4038917016667

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

