



Product description

| | |
|---------------------|--|
| Order number | 72.733.001 |
| Product description | Screw cap micro tube, working volume: 0.5 ml, conical base, with knurling, transparent, cap: natural, cap enclosed, 500 piece(s)/bag |

Product characteristics

| | |
|--------------------------|--------------------|
| Graduation | no |
| Cap | standard screw cap |
| Knurls | yes |
| Base shape | conical base |
| Volume of work | 0,5 ml |
| Centrifugation max (RCF) | 20000 x g |

Size

| | |
|-------------------|-------|
| Diameter | 8 mm |
| Length of product | 44 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

| | |
|-------------------|--|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |
| Closure material | Polypropylene (PP) |
| Colour of cap | natural |
| O-ring material | Ethylene propylene diene rubber (EPDM) |

Purity & certification

| | |
|---------------------------|-------------------------------|
| Satisfies the requirement | IATA |
| Product category | In vitrodiagnostic device, CE |
| CE certified | CE |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 5000 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 500 |
| Piece (s) / outer case | 5000 |
| Piece(s) / pallet | 180000 |
| Depth of case | 578 mm |
| Width of case | 368 mm |
| Height of case | 195 mm |
| Case volume | 0.0415 cbm |
| Weight of product | 0.0012 kg |
| Weight of case | 6.46 kg |
| EAN of inner box | 4038917119757 |
| EAN case | 4038917016773 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

