



Product description

Order number	72.733.201
Product description	Screw cap micro tube, for volume limitation to 100 µl, working volume: 0.5 ml, conical base, with knurling, transparent, without cap, 500 piece(s)/bag

Product characteristics

Graduation	no
Knurls	yes
Example of use	for volume limitation to 100 µl
Base shape	conical base
Volume of work	0,5 ml
Centrifugation max (RCF)	20000 x g

Size

Diameter	8 mm
Length of product	44 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

Screw cap micro tube, 0.5 ml

Page 2

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Satisfies the requirement	IATA
Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	500
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	5000
Piece(s) / pallet	210000
Depth of case	388 mm
Width of case	378 mm
Height of case	250 mm
Case volume	0.0367 cbm
Weight of product	0.0008 kg
Weight of case	4.48 kg
EAN of inner box	4038917366625
EAN case	4038917366618

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-03

This document was prepared by EDP support and is valid without signature.



SARSTEDT