



Product description

Order number	72.785
Product description	Screw cap micro tube, working volume: 0.5 ml, skirted conical base, transparent, without cap, 500 piece(s)/bag

Product characteristics

Graduation	no
Knurls	no
Base shape	skirted conical base
Volume of work	0,5 ml
Centrifugation max (RCF)	20000 x g

Size

Diameter	10.8 mm
Length of product	44 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

Screw cap micro tube, 0.5 ml

Page 2

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Satisfies the requirement	IATA
Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	5000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	5000
Piece(s) / pallet	180000
Depth of case	398 mm
Width of case	398 mm
Height of case	300 mm
Case volume	0.0475 cbm
Weight of product	0.0012 kg
Weight of case	6.63 kg
EAN of inner box	4038917084116
EAN case	4038917077675

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-20

This document was prepared by EDP support and is valid without signature.



SARSTEDT