



## Product Specification

**Tubo de rosca, 2,5 ml, (CxØ): 75 x 13 mm, fundo falso cônico, fundo do tubo arredondado, PP, sem tampa, 1.000 unid./pacote**

**Page 1**



### Product description

Order number	60.614.010
Descrição do artigo	Tubo de rosca, volume de trabalho: 2,5 ml, (CxØ): 75 x 13 mm, fundo falso cônico, fundo do tubo arredondado, transparente, material: PP, com escala, sem tampa, 1.000 unid./pacote, 2.000 unidades/caixa

### Product characteristics

Graduação	sim
Forma do fundo	fundo falso cônico, fundo do tubo arredondado
Volume de trabalho	2,5 ml

### Size

Diâmetro	13 mm
Altura do produto	75 mm
Comprimento do produto	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-04

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**Tubo de rosca, 2,5 ml, (CxØ): 75 x 13 mm, fundo falso cônico, fundo do tubo arredondado, PP, sem tampa, 1.000 unid./pacote**

**Page 2**

### Material & colours

Material do produto	Polipropileno (PP)
Cor do produto	transparente
Material do tubo	Polipropileno (PP)

### Purity & certification

Categoria do produto	Diagnóstico in vitro, CE
Certificação CE	CE
Em lotes	sim

### Packaging

Quantidade mínima do pedido	2000
Tipo de embalagem interior	pacote
Unidades por caixa interna	1000
Unidades por caixa	2000
Unidades por palete	48000
Profundidade da embalagem externa	458 mm
Largura da embalagem externa	408 mm
Altura da embalagem externa	260 mm
Volume da embalagem externa	0.0486 cbm
Peso do produto	0.0024 kg
Peso da embalagem externa	5.45 kg
EAN da caixa interna	4038917353182
EAN da embalagem externa	4038917353175

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-04

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**