



Product Specification

**S-Monovette® Soro CAT, 7,5 ml, tampa branca,
(CxØ): 92 x 15 mm, com etiqueta de papel**

Page 1



Product description

Order number	01.1601.014
Descrição do artigo	S-Monovette® Soro CAT, preparo: Ativador de coagulação, 7,5 ml, tampa de rosca com membrana, tampa branca, código de cor EU, (CxØ) sem tampa: 92 x 15 mm, com etiqueta de papel "Blutbank", 50 unid./caixa de papelão, estéril

Product characteristics

Tipo de retirada	venoso
Tipo de preparo	Ativador de coagulação
Código de cor	EU
Rótulo/impressão	com etiqueta de papel
Graduação	marcação nominal
Tampa	tampa de rosca com membrana
Tipo de fechamento	tampa de rosca
Centrifugação recomendada	2000 x g
Temperatura de centrifugação	20
Tempo de centrifugação	10 Min
Hermético	sim

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-02

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

**S-Monovette® Soro CAT, 7,5 ml, tampa branca,
(CxØ): 92 x 15 mm, com etiqueta de papel**

Page 2

Conexão tampa de rosca com membrana

Size

Volume nominal	7.5 ml
Diâmetro	15 mm
Comprimento com tampa	108 mm
Comprimento sem tampa	92 mm
Comprimento com tampa e barra	129 mm

Material & colours

Cor do produto	transparente
Material de fechamento	High Density Polyethylene (HD-PE)
Cor da tampa	branca
Material da membrana	borracha
Material do tubo	Polipropileno (PP)
Material do êmbolo	Poliétileno de alta densidade (HD-PE)
Material da haste do êmbolo	Poliestireno (PS)

Purity & certification

Corresponde às exigências	IATA, ADR
Categoria do produto	Diagnóstico in vitro, CE
Certificação CE	CE
Padrão de pureza	estéril
Esterilização	Irradiação de elétrons
Em lotes	sim

Packaging

Quantidade mínima do pedido	500
Tipo de embalagem interior	caixa de papelão
Unidades por caixa interna	50
Unidades por caixa	500
Unidades por palete	32000
Profundidade da caixa interna	189 mm
Largura da caixa interna	134 mm
Altura da caixa interna	82 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-02

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

**S-Monovette® Soro CAT, 7,5 ml, tampa branca,
(CxØ): 92 x 15 mm, com etiqueta de papel**

Page 3

Profundidade da embalagem externa	400 mm
Largura da embalagem externa	144 mm
Altura da embalagem externa	440 mm
Volume da embalagem externa	0.0253 cbm
Peso do produto	0.0084 kg
Peso da embalagem externa	5.04 kg
EAN da caixa interna	4038917131667
EAN da embalagem externa	4038917023849

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-02

This document was prepared by EDP support and is valid without signature.



SARSTEDT