



Product description

Order number	02.1063.001
Descrição do artigo	S-Monovette® Soro CAT, CAT, preparo: Ativador de coagulação, 9 ml, tampa de rosca com membrana, tampa branca, código de cor EU, (CxØ) sem tampa: 92 x 16 mm, com etiqueta de papel, rótulo/impressão: preta, 50 unid./caixa de papelão, estéril

Product characteristics

Tipo de retirada	venoso
Tipo de preparo	Ativador de coagulação
Código de cor	EU
Rótulo/impressão	com etiqueta de papel
Cor impressão/rótulo	preta
Graduação	marcação nominal
Tampa	tampa de rosca com membrana
Tipo de fechamento	tampa de rosca
Centrifugação recomendada	2000 x g
Temperatura de centrifugação	20
Tempo de centrifugação	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

**S-Monovette® Soro CAT, 9 ml, tampa branca,
(CxØ): 92 x 16 mm, com etiqueta de papel**

Page 2

Hermético sim
Conexão tampa de rosca com membrana

Size

Volume nominal 9 ml
Diâmetro externo do fecho 19 mm
Diâmetro 16 mm
Comprimento com tampa 108 mm
Comprimento sem tampa 92 mm
Comprimento com tampa e barra 129 mm

Material & colours

Cor do produto transparente
Material de fechamento High Density Polyethylene (HD-PE)
Cor da tampa branca
Material da membrana borracha
Material do tubo Polipropileno (PP)
Material do êmbolo Polietileno de alta densidade (HD-PE)
Material da haste do êmbolo Poliestireno (PS)

Purity & certification

Corresponde às exigências IATA, ADR
Categoria do produto Diagnóstico in vitro, CE
Certificação CE CE
Padrão de pureza estéril
Esterilização Irradiação de elétrons
Em lotes sim

Packaging

Quantidade mínima do pedido 500
Tipo de embalagem interior caixa de papelão
Unidades por caixa interna 50
Unidades por caixa 500
Unidades por palete 24000
Profundidade da caixa interna 204 mm
Largura da caixa interna 134 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-27

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
S-Monovette® Soro CAT, 9 ml, tampa branca,
(CxØ): 92 x 16 mm, com etiqueta de papel

Page 3

Altura da caixa interna	88 mm
Profundidade da embalagem externa	648 mm
Largura da embalagem externa	303 mm
Altura da embalagem externa	150 mm
Volume da embalagem externa	0.0295 cbm
Peso do produto	0.0088 kg
Peso da embalagem externa	5.35 kg
EAN da caixa interna	4038917124959
EAN da embalagem externa	4038917019569

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-27

This document was prepared by EDP support and is valid without signature.



SARSTEDT