



Product Specification

Recipiente de envio refrigerado para capilares de gases sanguíneos, S-Monovette® até 105 x Ø 18 mm, transparente, comprimento: 50 mm

Page 1



Product description

Order number	95.997
Descrição do artigo	Recipiente de envio refrigerado para capilares de gases sanguíneos, S-Monovette® até 105 x Ø 18 mm, transparente, comprimento: 50 mm, Transporte refrigerado/transporte aquecido, 1 unid./caixa de papelão

Product characteristics

Campo de aplicação Indicado para	Transporte refrigerado/transporte aquecido capilares de gases sanguíneos, S-Monovette® até 105 x Ø 18 mm
Compressa de absorção	não

Size

Largura do produto	50 mm
Altura do produto	140 mm
Comprimento do produto	50 mm
Comprimento útil	118 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-20

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

Recipiente de envio refrigerado para capilares de gases sanguíneos, S-Monovette® até 105 x Ø 18 mm, transparente, comprimento: 50 mm

Page 2

Material & colours

Material do produto	Polietileno de alta densidade (HD-PE)
Cor do produto	transparente
Cor da tampa	transparente

Purity & certification

Categoria do produto	Não é um dispositivo médico DIV
Em lotes	sim

Packaging

Quantidade mínima do pedido	1
Tipo de embalagem interior	caixa de papelão
Unidades por caixa interna	1
Unidades por caixa	66
Unidades por palete	2508
Profundidade da embalagem externa	568 mm
Largura da embalagem externa	433 mm
Altura da embalagem externa	120 mm
Volume da embalagem externa	0.0295 cbm
Peso do produto	0.0446 kg
Peso da embalagem externa	4.59 kg
EAN da caixa interna	4038917084918
EAN da embalagem externa	4038917362382

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-20

This document was prepared by EDP support and is valid without signature.



SARSTEDT