



## Product Specification

**Embalagem de transporte, LabMailer Small, com compressa de absorção, comprimento: 82 mm, Ø da abertura: 78 mm, tampa montada**

**Page 1**



### Product description

Order number	78.9925.418
Descrição do artigo	Embalagem de transporte, LabMailer Small, com compressa de absorção, (CxL): 82 x 82 mm, Ø da abertura: 78 mm, tampa montada, tampa: vermelha

### Product characteristics

Campo de aplicação	Transporte de amostras para diagnóstico
Compressa de absorção	sim

### Size

Diâmetro	82 mm
Diâmetro da abertura	78 mm
Largura do produto	82 mm
Comprimento do produto	82 mm
Comprimento com tampa	179 mm
Comprimento útil	153 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-26

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**Embalagem de transporte, LabMailer Small, com compressa de absorção, comprimento: 82 mm, Ø da abertura: 78 mm, tampa montada**

**Page 2**

### Material & colours

Material do produto	papelão
Cor da tampa	vermelha

### Purity & certification

Isento de látex	sim
Isenção de pirogênio	não
Ficha de dados de segurança	não

### Packaging

Quantidade mínima do pedido	50
Unidades por caixa interna	50
Unidades por caixa	50
Unidades por palete	1600
Volume da embalagem externa	0.0692 cbm
Peso da embalagem externa	11 kg
EAN da caixa interna	4038917155526
EAN da embalagem externa	4038917039437

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-26

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**