



Product Specification

UriSet 24, sem preparação, 3 l, com Monovette® para urina, marrom, com proteção contra a luz, graduado

Page 1



Product description

| | |
|---------------------|--|
| Order number | 77.578.010 |
| Descrição do artigo | UriSet 24, Conjunto completo para coleta ideal de urina de 24 horas, volume de trabalho máx.: 3 l, Ø da abertura: 83,5 mm, sem preparação, marrom, com proteção contra a luz, graduado, material: PE, tampa: amarela, com Monovette® para urina, 30 unid./caixa de papelão |

Product characteristics

| | |
|------------------------|--|
| Tipo de preparo | sem preparação |
| Graduação | sim |
| Tipo de fechamento | tampa de rosca |
| Proteção contra a luz | sim |
| Versão | com Monovette® para urina |
| Campo de aplicação | Conjunto completo para coleta ideal de urina de 24 horas |
| Janela de visualização | sim |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2024-03-29

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

**UriSet 24, sem preparação, 3 l, com Monovette®
para urina, marrom, com proteção contra a luz,
graduado**

Page 2

Size

| | |
|------------------------|-------------|
| Volume nominal | 3.000.00 ml |
| Volume máximo | 3 l |
| Diâmetro da abertura | 83,5 mm |
| Largura do produto | 130 mm |
| Comprimento do produto | 125 mm |
| Comprimento com tampa | 246 mm |
| Comprimento sem tampa | 243 mm |

Material & colours

| | |
|------------------------|--------------------|
| Material do produto | Polietileno (PE) |
| Cor do produto | marrom |
| Material de fechamento | Polypropylene (PP) |
| Cor da tampa | amarela |

Purity & certification

| | |
|----------------------|--------------------------|
| Categoria do produto | Diagnóstico in vitro, CE |
| Certificação CE | CE |
| Em lotes | sim |

Packaging

| | |
|--------------------------------------|------------------|
| Quantidade mínima do pedido | 30 |
| Tipo de embalagem interior | caixa de papelão |
| Unidades por caixa interna | 30 |
| Unidades por caixa | 30 |
| Unidades por palete | 270 |
| Profundidade da embalagem externa | 643 mm |
| Largura da embalagem externa | 403 mm |
| Altura da embalagem externa | 515 mm |
| Volume da embalagem externa | 0.1335 cbm |
| Peso do produto | 0.1846 kg |
| Peso da embalagem externa | 6.34 kg |
| EAN da caixa interna | 4038917113861 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2024-03-29

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

**UriSet 24, sem preparação, 3 l, com Monovette®
para urina, marrom, com proteção contra a luz,
graduado**

Page 3

EAN da embalagem externa 4038917080019

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2024-03-29

This document was prepared by EDP support and is valid without signature.



SARSTEDT