



Product Specification

Rotor de placas, para 2 Tubos de até 35 mm de Ø e 6 tubos de até 20 mm de Ø e 6 tubos de até 12,5 mm de Ø, para SARMIX® M 2000

Page 1



Product description

Order number	92.180.611
Descrição do artigo	Rotor de placas, para 2 Tubos de até 35 mm de Ø e 6 tubos de até 20 mm de Ø e 6 tubos de até 12,5 mm de Ø, para SARMIX® M 2000, 1 unid./caixa de papelão

Product characteristics

Indicado para	SARMIX® M 2000
Capacidade de armazenamento	2

Purity & certification

Categoria do produto	Não é um dispositivo médico DIV
Em lotes	não

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

Rotor de placas, para 2 Tubos de até 35 mm de Ø e 6 tubos de até 20 mm de Ø e 6 tubos de até 12,5 mm de Ø, para SARMIX® M 2000

Page 2

Packaging

Quantidade mínima do pedido	1
Tipo de embalagem interior	caixa de papelão
Unidades por caixa interna	1
Unidades por caixa	10
Unidades por palete	10
Volume da embalagem externa	0.0122 cbm
Peso do produto	0.2009 kg
Peso da embalagem externa	2.21 kg
EAN da caixa interna	4038917101998
EAN da embalagem externa	4038917075756

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-26

This document was prepared by EDP support and is valid without signature.



SARSTEDT