



Product description

Order number	92.181.611
Descrição do artigo	Rotor, para 22 recipientes

Product characteristics

Capacidade de armazenamento	22
-----------------------------	----

Purity & certification

Categoria do produto	não regulamentado
Certificação CE	sem marca CE
Em lotes	sim

Packaging

Quantidade mínima do pedido	1
Unidades por caixa interna	1
Unidades por caixa	1
Unidades por palete	40
Volume da embalagem externa	0.0052 cbm
Peso do produto	1.3595 kg

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Rotor, para 22 recipientes

Page 2

Peso da embalagem externa 1.36 kg
EAN da caixa interna 4038917156257
EAN da embalagem externa 4038917156240

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2026-03-07

This document was prepared by EDP support and is valid without signature.



SARSTEDT