



**Product Specification**  
**Tubo CryoPure, 1,2 ml, tampa de rosca**  
**QuickSeal, cores variadas**

**Page 1**



**Product description**

Order number	72.377.992
Descrição do artigo	Tubo CryoPure, 1,2 ml, tubo: PP, tampa de rosca QuickSeal, tampa montada, HD-PE, cores variadas, rosca externa, volume nominal: 1 ml, Cryo Performance Tested, 50 unid./pacote

**Product characteristics**

Rótulo/impressão	com impressão
Cor impressão/rótulo	branca
Graduação	sim
Tampa	tampa de rosca QuickSeal
Tipo de fechamento	tampa de rosca
Rosca	rosca externa
Volume de trabalho	1 ml
Tipo de recipiente	Tubo

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-27

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



### Size

Volume nominal	1 ml
Volume máximo	1.2 ml
Diâmetro	12.65 mm
Comprimento do produto	36.7 mm
Comprimento sem tampa	29.8 mm

### Material & colours

Material de fechamento	High Density Polyethylene (HD-PE)
Cor da tampa	cores variadas
Material do tubo	Polipropileno (PP)

### Purity & certification

Categoria do produto	Diagnóstico in vitro, CE
Certificação CE	CE
Padrão de pureza	Cryo Performance Tested
Esterilização	Irradiação de elétrons
Isenção de pirogênio	sim
Não citotóxico	sim
Em lotes	sim

### Packaging

Quantidade mínima do pedido	500
Tipo de embalagem interior	pacote
Unidades por caixa interna	500
Unidades por caixa	2000
Unidades por palete	104000
Profundidade da caixa interna	293 mm
Largura da caixa interna	183 mm
Altura da caixa interna	118 mm
Profundidade da embalagem externa	598 mm
Largura da embalagem externa	378 mm
Altura da embalagem externa	135 mm
Volume da embalagem externa	0.0305 cbm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Tubo CryoPure, 1,2 ml, tampa de rosca**  
**QuickSeal, cores variadas**

**Page 3**

Peso do produto	0.0015 kg
Peso da embalagem externa	4.25 kg
EAN da caixa interna	4038917228428
EAN da embalagem externa	4038917228411

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-27

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**