



## Product Specification

**Copo, Tiosulfato de sódio, 250 ml, (CxØ): 78 x 70 mm, graduado, PS, transparente, com etiqueta de papel**

Page 1



### Product description

Order number	51.9923.821
Descrição do artigo	Copo, Tiosulfato de sódio, volume de trabalho máx.: 250 ml, (CxØ): 78 x 70 mm, transparente, tampa: amarela, graduado, material: PS, tampa de rosca, tampa montada, com etiqueta de papel, rótulo/impressão: preta, estéril, 210 unid./pacote

### Product characteristics

Tipo de preparo	Solução de tiosulfato de sódio
Rótulo/impressão	com etiqueta de papel
Cor impressão/rótulo	preta
Graduação	sim
Tampa	tampa de rosca padrão
Tipo de fechamento	tampa de rosca
Campo de aplicação	Recipiente para amostras de água
Forma do fundo	fundo plano
Volume de trabalho	250 ml

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-20

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**Copo, Tiosulfato de sódio, 250 ml, (CxØ): 78 x 70 mm, graduado, PS, transparente, com etiqueta de papel**

Page 2

### Size

Volume nominal	250 ml
Diâmetro	70 mm
Altura do produto	78 mm
Comprimento do produto	78 mm

### Material & colours

Material do produto	Poliestireno (PS)
Cor do produto	transparente
Material de fechamento	High Density Polyethylene (HD-PE)
Cor da tampa	amarela
Material do tubo	Poliestireno (PS)

### Purity & certification

Corresponde às exigências	hermético a 9"hg/10 min
Isento de látex	sim
Padrão de pureza	estéril
Isenção de pirogênio	não
Ficha de dados de segurança	não

### Packaging

Quantidade mínima do pedido	210
Tipo de embalagem interior	pacote
Unidades por caixa interna	240
Unidades por caixa	210
Unidades por palete	5040
Profundidade da embalagem externa	534 mm
Largura da embalagem externa	384 mm
Altura da embalagem externa	519 mm
Volume da embalagem externa	0.1025 cbm
Peso da embalagem externa	10.15 kg
EAN da caixa interna	4038917154567
EAN da embalagem externa	4038917235068

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-20

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**