



Product description

Order number	70.1186.210
Descrição do artigo	Ponta de filtro, volume de trabalho: 1.250 µl, transparente, anéis de nível de enchimento, Biosphere® , adequado para SARSTEDT Sarpette® M, Eppendorf, Gilson e Brand e desenhos idênticos, 96 unid./caixa

Product characteristics

Graduação	sim
Indicado para	SARSTEDT Sarpette® M, Eppendorf, Gilson e Brand, bem como modelos idênticos
Volume de trabalho	1.250 µl
Filtro	sim
Cor de codificação	azul

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Ponta de filtro, 1.250 µl, transparente,
Biosphere® , 96 unid./caixa

Page 2

Size

Volume máximo	1250 µl
Comprimento do produto	103 mm

Material & colours

Material do produto	Polipropileno (PP)
Cor do produto	transparente
Material do filtro	Polietileno (PE)

Purity & certification

Categoria do produto	Não é um dispositivo médico DIV
Padrão de pureza	Biosphere®
Esterilização	Radiação gama
Isenção de pirogênio	sim
Em lotes	sim

Packaging

Quantidade mínima do pedido	768
Tipo de embalagem interior	caixa
Unidades por caixa interna	96
Unidades por caixa	768
Unidades por palete	90000
Profundidade da embalagem externa	270 mm
Largura da embalagem externa	121 mm
Altura da embalagem externa	413 mm
Volume da embalagem externa	0.0135 cbm
Peso do produto	0.0022 kg
Peso da embalagem externa	2.08 kg
EAN da caixa interna	4038917126496
EAN da embalagem externa	4038917080354

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2024-03-29

This document was prepared by EDP support and is valid without signature.



SARSTEDT