



Product description

Order number	70.1193.151
Descrição do artigo	Ponta de filtro, volume de trabalho: 1.000 µl, preta, adequado para aplicação em dispositivos SARSTEDT IVARO AP, 480 unid./blister

Product characteristics

Indicado para	aplicação em dispositivos SARSTEDT IVARO AP
Volume de trabalho	1.000 µl
Filtro	sim

Size

Volume máximo	1000 µl
---------------	---------

Material & colours

Material do produto	Polipropileno (PP)
Cor do produto	preta

Purity & certification

Categoria do produto	Não é um dispositivo médico DIV
Em lotes	sim

Packaging

Quantidade mínima do pedido	3840
Tipo de embalagem interior	blister
Unidades por caixa interna	480
Unidades por caixa	3840
Unidades por palete	61440
Profundidade da embalagem externa	500 mm
Largura da embalagem externa	335 mm
Altura da embalagem externa	200 mm
Volume da embalagem externa	0.0754 cbm
Peso do produto	0.002 kg
Peso da embalagem externa	7.95 kg
EAN da caixa interna	4038917913447

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Ponta de filtro, 1.000 µl, preta, 480 unid./blister

Page 2

EAN da embalagem externa 4038917913454

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-16

This document was prepared by EDP support and is valid without signature.



SARSTEDT