




## Product Specification

**Extremidade de pipeta, 1.000 µl, transparente,  
PCR Performance Tested, Low Retention, 96  
unid./caixa**

**Page 1**



### Product description

Order number	70.3050.210
Descrição do artigo	Extremidade de pipeta, volume de trabalho: 1.000 µl, transparente, anéis de nível de enchimento, PCR Performance Tested, Low Retention, adequado para SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit e Brand, bem como modelos idênticos, 96 unid./caixa

### Product characteristics

Gradação	sim
Indicado para	SARSTEDT Sarpette® M, Eppendorf, Gilson, Finnpipette, Biohit e Brand, bem como modelos idênticos
Volume de trabalho	1.000 µl
Filtro	não
Cor de codificação	azul
Autoclavável (máx. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2026-04-04

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**Extremidade de pipeta, 1.000 µl, transparente,  
PCR Performance Tested, Low Retention, 96  
unid./caixa**

**Page 2**

### Size

Volume máximo	1000 µl
Diâmetro	9 mm
Comprimento do produto	72 mm

### Material & colours

Material do produto	Polipropileno (PP)
Cor do produto	transparente

### Purity & certification

Categoria do produto	Diagnóstico in vitro (IVD)
Certificado de conformidade	sim
Certificação CE	CE - Declaração do fabricante
Padrão de pureza	PCR Performance Tested
Em lotes	sim

### Packaging

Quantidade mínima do pedido	1920
Tipo de embalagem interior	caixa
Unidades por caixa interna	480
Unidades por caixa	1920
Unidades por palete	61440
Profundidade da caixa interna	530 mm
Largura da caixa interna	111 mm
Altura da caixa interna	148 mm
Profundidade da embalagem externa	547 mm
Largura da embalagem externa	325 mm
Altura da embalagem externa	245 mm
Volume da embalagem externa	0.0436 cbm
Peso do produto	0.0018 kg
Peso da embalagem externa	4.95 kg
EAN da caixa interna	4038917424721
EAN da embalagem externa	4038917424714

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2026-04-04

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**