



Product description

Order number	86.1173.300
Descrição do artigo	Pipeta de transferência, volume nominal: 3,5 ml, (CxL): 155 x 12,5 mm, material: PEBD, transparente, tamanho da gota: 30 µl - 45 µl, 500 unid./caixa de papelão

Product characteristics

Graduação	não
-----------	-----

Size

Volume nominal	3.5 ml
Largura do produto	12.5 mm
Comprimento do produto	155 mm
Tamanho da gota	30 µl - 45 µl

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Pipeta de transferência, 3,5 ml, (CxL): 155 x 12,5 mm, PEBD, transparente

Page 2

Material & colours

Material do produto Polietileno de baixa densidade (PEBD)
Cor do produto transparente

Purity & certification

Categoria do produto Não é um dispositivo médico | DIV
Em lotes sim

Packaging

Quantidade mínima do pedido 500
Tipo de embalagem interior caixa de papelão
Unidades por caixa interna 500
Unidades por caixa 5000
Unidades por palete 56000
Volume da embalagem externa 0.0994 cbm
Peso do produto 0.0012 kg
Peso da embalagem externa 6.89 kg
EAN da caixa interna 4038917378475
EAN da embalagem externa 4038917378468

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-10-06

This document was prepared by EDP support and is valid without signature.



SARSTEDT