



Product description

Order number	86.1175.300
Descrição do artigo	Pipeta de transferência, volume nominal: 6 ml, (CxL): 146 x 15 mm, material: PEBD, transparente, tamanho da gota: 20 µl - 30 µl

Product characteristics

Graduação	não
-----------	-----

Size

Volume nominal	6 ml
Largura do produto	15 mm
Comprimento do produto	146 mm
Tamanho da gota	20 µl - 30 µl

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Pipeta de transferência, 6 ml, (CxL): 146 x 15 mm,
PEBD, transparente

Page 2

Material & colours

Material do produto Polietileno de baixa densidade (PEBD)
Cor do produto transparente

Purity & certification

Categoria do produto Não é um dispositivo médico | DIV
Em lotes sim

Packaging

Quantidade mínima do pedido 400
Unidades por caixa interna 400
Unidades por caixa 4000
Unidades por palete 44800
Volume da embalagem 0.0994 cbm
externa
Peso do produto 0.0018 kg
Peso da embalagem externa 7.95 kg
EAN da caixa interna 4038917378536
EAN da embalagem externa 4038917378529

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-26

This document was prepared by EDP support and is valid without signature.



SARSTEDT