




Product description

Order number	95.2092
Descrição do artigo	Filme, livre de Dnase/Rnase, transparente, adequado para qPCR, perfurável, campo de aplicação: -80 - 120 °C

Product characteristics

Tipo de fechamento	filme
Proteção contra a luz	não
Indicada para qPCR	sim
Temperatura mínima	-80 °C
Temperatura máxima	120 °C

Size

Largura do produto	82 mm
Comprimento do produto	125 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Filme, livre de Dnase/Rnase, transparente,
perfurável**

Page 2

Material & colours

Cor do produto transparente

Purity & certification

Categoria do produto Não é um dispositivo médico | DIV
Padrão de pureza livre de Dnase/Rnase
Em lotes sim

Packaging

Quantidade mínima do pedido 100
Unidades por caixa interna 100
Unidades por caixa 1200
Unidades por palete 57600
Volume da embalagem 0.0142 cbm
externa
Peso do produto 0.0031 kg
Peso da embalagem externa 4.36 kg
EAN da caixa interna 4038917335676
EAN da embalagem externa 4038917335669

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-03

This document was prepared by EDP support and is valid without signature.



SARSTEDT