



### **Product description**

Order number	95.64.259
Descrição do artigo	Rack, material: poliestireno expandido, branca, dimensão da grade: 5 x 4, para 20 recipientes, 8 unid./caixa de papelão

### **Product characteristics**

Capacidade de armazenamento	20
Dimensão da grade	5 x 4

### **Material & colours**

Material do produto	poliestireno expandido
Cor do produto	branca

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Rack, poliestireno expandido, dimensão da grade: 5 x 4**

**Page 2**

**Purity & certification**

Categoria do produto Não é um dispositivo médico | DIV  
Em lotes sim

**Packaging**

Quantidade mínima do pedido 56  
Tipo de embalagem interior caixa de papelão  
Unidades por caixa interna 8  
Unidades por caixa 56  
Unidades por palete 1792  
Volume da embalagem externa 0.0559 cbm  
Peso do produto 0.0238 kg  
Peso da embalagem externa 1.36 kg  
EAN da caixa interna 4038917386883  
EAN da embalagem externa 4038917386876

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-10-06

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**