




Product description

Order number	55.526.002
Descrição do artigo	Tubo, volume de trabalho: 5 ml, (CxØ): 75 x 12 mm, material: PP, fundo redondo, marrom, tampa de pressão, 500 unid./embalagem empilhável

Product characteristics

Graduação	não
Tipo de fechamento	tampa de pressão
Forma do fundo	fundo redondo
Volume de trabalho	5 ml
Centrifugação máxima (RCF)	5000 x g
Autoclavável (máx. °C)	121 °C

Size

Diâmetro	12 mm
Altura do produto	75 mm
Comprimento do produto	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Tubo, 5 ml, (CxØ): 75 x 12 mm, PP

Page 2

Material & colours

Material do produto	Polipropileno (PP)
Cor do produto	marrom

Purity & certification

Categoria do produto	Não é um dispositivo médico DIV
Em lotes	sim

Packaging

Quantidade mínima do pedido	4000
Tipo de embalagem interior	embalagem empilhável
Unidades por caixa interna	500
Unidades por caixa	4000
Unidades por palete	120000
Volume da embalagem externa	0.046 cbm
Peso do produto	0.0021 kg
Peso da embalagem externa	10.48 kg
EAN da caixa interna	4038917109284
EAN da embalagem externa	4038917031974

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT