



Product Specification

Tubo de rosca, 5 ml, (CxØ): 75 x 13 mm, fundo redondo, PP, tampa incluída, 100 unid./pacote

Page 1



Product description

Order number	62.504.039
Descrição do artigo	Tubo de rosca, volume de trabalho: 5 ml, (CxØ): 75 x 13 mm, fundo redondo, marrom, material: PP, com impressão, rótulo/impressão: preta, com escala, tampa incluída, marrom, 100 unid./pacote, 1.000 unidades/caixa

Product characteristics

Rótulo/impressão	com impressão
Cor impressão/rótulo	preta
Graduação	sim
Tipo de fechamento	tampa de rosca
Proteção contra a luz	sim
Forma do fundo	fundo redondo
Volume de trabalho	5 ml

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

Tubo de rosca, 5 ml, (CxØ): 75 x 13 mm, fundo redondo, PP, tampa incluída, 100 unid./pacote

Page 2

Size

Diâmetro	13 mm
Comprimento do produto	75 mm

Material & colours

Material do produto	Polipropileno (PP)
Cor do produto	marrom
Material de fechamento	High Density Polyethylene (HD-PE)
Cor da tampa	marrom
Material do tubo	Polipropileno (PP)

Purity & certification

Categoria do produto	Diagnóstico in vitro, CE
Certificação CE	CE
Em lotes	sim

Packaging

Quantidade mínima do pedido	50000
Tipo de embalagem interior	pacote
Unidades por caixa interna	100
Unidades por caixa	1000
Unidades por palete	40000
Profundidade da embalagem externa	458 mm
Largura da embalagem externa	408 mm
Altura da embalagem externa	205 mm
Volume da embalagem externa	0.0383 cbm
Peso do produto	0.0039 kg
Peso da embalagem externa	4.39 kg
EAN da caixa interna	4038917413008
EAN da embalagem externa	4038917412995

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT