



## Product Specification

Tubo de rosca, 5 ml, (CxØ): 92 x 15,3 mm, fundo falso cônico, fundo do tubo plano, PP, tampa montada, 100 unid./pacote

Page 1



### Product description

Order number	62.611
Descrição do artigo	Tubo de rosca, volume de trabalho: 5 ml, (CxØ): 92 x 15,3 mm, fundo falso cônico, fundo do tubo plano, transparente, material: PP, com impressão, rótulo/impressão: branca, com escala, tampa montada, 100 unid./pacote, 1.000 unidades/caixa

### Product characteristics

Rótulo/impressão	com impressão
Cor impressão/rótulo	branca
Graduação	sim
Forma do fundo	fundo falso cônico, fundo do tubo plano
Volume de trabalho	5 ml

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**Tubo de rosca, 5 ml, (CxØ): 92 x 15,3 mm, fundo falso cônico, fundo do tubo plano, PP, tampa montada, 100 unid./pacote**

**Page 2**

### Size

Diâmetro	15.3 mm
Comprimento do produto	92 mm
Comprimento sem tampa	92 mm

### Material & colours

Material do produto	Polipropileno (PP)
Cor do produto	transparente

### Purity & certification

Corresponde às exigências	IATA, ADR
Categoria do produto	Diagnóstico in vitro, CE
Certificação CE	CE
Em lotes	sim

### Packaging

Quantidade mínima do pedido	1000
Tipo de embalagem interior	pacote
Unidades por caixa interna	100
Unidades por caixa	1000
Unidades por palete	36000
Profundidade da embalagem externa	398 mm
Largura da embalagem externa	398 mm
Altura da embalagem externa	300 mm
Volume da embalagem externa	0.0475 cbm
Peso do produto	0.005 kg
Peso da embalagem externa	5.69 kg
EAN da caixa interna	4038917082877
EAN da embalagem externa	4038917001892

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**