



Product Specification

Copo multiuso, 250 ml, (CxØ): 78 x 70 mm, graduado, PP, transparente, com rótulo de segurança

Page 1



Product description

| | |
|---------------------|--|
| Order number | 75.9922.518 |
| Descrição do artigo | Copo multiuso, volume de trabalho máx.: 250 ml, (CxØ): 78 x 70 mm, transparente, tampa: amarela, graduado, material: PP, tampa de rosca, tampa montada, com rótulo de segurança, estéril, 240 unid./pacote |

Product characteristics

| | |
|--------------------|-------------------------|
| Rótulo/impressão | com rótulo de segurança |
| Graduação | sim |
| Tipo de fechamento | tampa de rosca |
| Forma do fundo | fundo plano |
| Volume de trabalho | 250 ml |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

Copo multiuso, 250 ml, (CxØ): 78 x 70 mm, graduado, PP, transparente, com rótulo de segurança

Page 2

Size

| | |
|------------------------|--------|
| Volume nominal | 250 ml |
| Diâmetro | 70 mm |
| Altura do produto | 78 mm |
| Comprimento do produto | 78 mm |

Material & colours

| | |
|------------------------|-----------------------------------|
| Material do produto | Polipropileno (PP) |
| Cor do produto | transparente |
| Material de fechamento | High Density Polyethylene (HD-PE) |
| Cor da tampa | amarela |

Purity & certification

| | |
|----------------------|-----------------------------------|
| Categoria do produto | Não é um dispositivo médico DIV |
| Padrão de pureza | estéril |
| Esterilização | Radiação gama |
| Em lotes | sim |

Packaging

| | |
|-----------------------------------|---------------|
| Quantidade mínima do pedido | 240 |
| Tipo de embalagem interior | pacote |
| Unidades por caixa interna | 240 |
| Unidades por caixa | 240 |
| Unidades por palete | 1440 |
| Profundidade da embalagem externa | 646 mm |
| Largura da embalagem externa | 450 mm |
| Altura da embalagem externa | 498 mm |
| Volume da embalagem externa | 0.1448 cbm |
| Peso do produto | 0.0404 kg |
| Peso da embalagem externa | 10.93 kg |
| EAN da caixa interna | 4038917250955 |
| EAN da embalagem externa | 4038917350914 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT