




**Product Specification**  
**Blood gas capillaries, calcium-balanced lithium  
heparin, 175 µl**

**Page 1**



**Product description**

Order number	19.930.175
Product description	Blood gas capillaries, preparation: calcium-balanced lithium heparin, nominal volume: 175 µl, (LxØ): 100 x 2.3 mm, 200 piece(s)/container

**Product characteristics**

Type of collection	capillary
Type of preparation	calcium-balanced lithium heparin
Preparation concentration	70 IU/ml blood
Application	blood gas analysis

**Size**

Sample volume	175 µl
Diameter	2.3 mm
Length of product	100 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

Product material	Polyethylene therephthalate (PET)
Colour of product	transparent

**Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Latex-free	yes
Batched	yes

**Packaging**

Minimum order qty.	200
Type of smallest subpackaging	container
Piece(s) / inner box	200
Piece (s) / outer case	12800
Piece(s) / pallet	921600
Depth of case	333 mm
Width of case	303 mm
Height of case	190 mm
Case volume	0.0192 cbm
Weight of product	0.0003 kg
Weight of case	5.5 kg
EAN of inner box	4038917234344
EAN case	4038917234337

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

