



Product description

| | |
|---------------------|---|
| Order number | 80.9924.014 |
| Product description | Faeces tube, with blade, white, screw cap, cap: brown, (LxØ) without cap: 54 x 28 mm, with paper label, 50 piece(s)/bag |

Product characteristics

| | |
|-----------------------|------------------|
| Label/ Print | with paper label |
| Colour of print/label | black |
| Closure type | screw cap |
| Design | blade |

Size

| | |
|----------------------|-------|
| Diameter | 28 mm |
| Length of product | 55 mm |
| Length excluding cap | 54 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

Faeces tube, with blade, screw cap, (LxØ): 54 x 28 mm, white

Page 2

Material & colours

| | |
|-------------------|-----------------------------------|
| Colour of product | white |
| Closure material | High Density Polyethylene (HD-PE) |
| Colour of cap | brown |
| Tube material | Polypropylene (PP) |

Purity & certification

| | |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified | CE |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 500 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 50 |
| Piece (s) / outer case | 500 |
| Piece(s) / pallet | 14000 |
| Depth of case | 458 mm |
| Width of case | 373 mm |
| Height of case | 238 mm |
| Case volume | 0.0407 cbm |
| Weight of product | 0.0084 kg |
| Weight of case | 4.9 kg |
| EAN of inner box | 4038917154772 |
| EAN case | 4038917039598 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-06

This document was prepared by EDP support and is valid without signature.



SARSTEDT