



Product Specification
**Multivette® 600 Serum CAT, 600 µl, cap white,
screw cap**

Page 1



Product description

| | |
|---------------------|---|
| Order number | 15.1670 |
| Product description | Multivette® 600 Serum CAT, venous and capillary blood samples, preparation: Clotting activator, nominal volume: 600 µl, (LxØ) with cap: 47.9 x 8 mm, with print, cap white, screw cap, flat base, 100 piece(s)/case |

Product characteristics

| | |
|----------------------------|------------------------------------|
| Type of collection | Venous and capillary |
| Type of preparation | Clotting activator |
| Colour code | EU |
| Label/ Print | with print |
| Colour of print/label | white |
| Cap | standard screw cap |
| Closure type | screw cap |
| Application | venous and capillary blood samples |
| Base shape | flat base |
| Recommended centrifugation | 10000 U/min |
| Centrifugation time | 5 Min |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Size

| | |
|----------------------|---------|
| Sample volume | 600 µl |
| Diameter | 8 mm |
| Length including cap | 47.9 mm |
| Length excluding cap | 44.4 mm |

Material & colours

| | |
|-------------------|-----------------------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |
| Closure material | High Density Polyethylene (HD-PE) |
| Colour of cap | white |

Purity & certification

| | |
|------------------|--------------------------------|
| Product category | In vitro diagnostic device, CE |
| CE certified | CE |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1000 |
| Type of smallest subpackaging | case |
| Piece(s) / inner box | 100 |
| Piece (s) / outer case | 1000 |
| Piece(s) / pallet | 63000 |
| Depth of box | 201 mm |
| Width of box | 81 mm |
| Height of box | 81 mm |
| Depth of case | 423 mm |
| Width of case | 213 mm |
| Height of case | 185 mm |
| Case volume | 0.0167 cbm |
| Weight of product | 0.0021 kg |
| Weight of case | 2.76 kg |
| EAN of inner box | 4038917087384 |
| EAN case | 4038917005814 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

