



Product description

| | |
|---------------------|--|
| Order number | 16.443.100 |
| Product description | Microvette® CB 300 Lithium heparin LH, Veterinary medicine, preparation: Lithium heparin, nominal volume: 300 µl, (LxØ) with cap: 46.5 x 10.8 mm, cap green, push cap, flat base, inner receptacle: trumpet-shaped, 100 piece(s)/bag |

Product characteristics

| | |
|----------------------------|---------------------|
| Type of collection | capillary |
| Type of preparation | Lithium heparin |
| Colour code | ISO |
| Closure type | push cap |
| Application | Veterinary medicine |
| Base shape | flat base |
| Recommended centrifugation | 2000 U/min |
| Centrifugation time | 5 Min |
| Shape of inner tube | trumpet-shaped |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Microvette® CB 300 Lithium heparin LH, 300 µl,
cap green, push cap, flat base**

Page 2

Size

| | |
|----------------------|---------|
| Sample volume | 300 µl |
| Diameter | 10.8 mm |
| Length including cap | 46.5 mm |
| Length excluding cap | 43.7 mm |

Material & colours

| | |
|-------------------|----------------------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |
| Closure material | Low Density Polyethylene (LD-PE) |
| Colour of cap | green |

Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1000 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 100 |
| Piece (s) / outer case | 1000 |
| Piece(s) / pallet | 80000 |
| Depth of case | 378 mm |
| Width of case | 288 mm |
| Height of case | 180 mm |
| Case volume | 0.0196 cbm |
| Weight of product | 0.002 kg |
| Weight of case | 2.35 kg |
| EAN of inner box | 4038917268349 |
| EAN case | 4038917282765 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

