



Product description

| | |
|---------------------|---|
| Order number | 16.443 |
| Product description | Microvette® CB 300 Lithium heparin LH, capillary blood collection, preparation: lithium heparin LH, nominal volume: 300 µl, (LxØ) with cap: 46.5 x 10.8 mm, cap orange, push cap, flat base, inner receptacle: trumpet-shaped, 100 piece(s)/bag |

Product characteristics

| | |
|----------------------------|----------------------------|
| Type of collection | capillary |
| Type of preparation | lithium heparin LH |
| Colour code | EU |
| Cap | Softstopfen mit Schlaufe |
| Closure type | push cap |
| Application | capillary blood collection |
| Base shape | flat base |
| Recommended centrifugation | 2000 U/min |
| Centrifugation time | 5 Min |
| Shape of inner tube | trumpet-shaped |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Size

| | |
|----------------------|---------|
| Sample volume | 300 µl |
| Diameter | 10.8 mm |
| Length including cap | 46.5 mm |
| Length excluding cap | 43.7 mm |

Material & colours

| | |
|-------------------|----------------------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |
| Closure material | Low Density Polyethylene (LD-PE) |
| Colour of cap | orange |

Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 1000 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 100 |
| Piece (s) / outer case | 1000 |
| Piece(s) / pallet | 80000 |
| Depth of case | 378 mm |
| Width of case | 288 mm |
| Height of case | 180 mm |
| Case volume | 0.0196 cbm |
| Weight of product | 0.002 kg |
| Weight of case | 2.35 kg |
| EAN of inner box | 4038917268295 |
| EAN case | 4038917282727 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

