



Product Specification

**Minivette® POCT Lithium heparin LH, 200 µl,
plunger green, colour code ISO, 150 piece(s)/bag**

Page 1



Product description

Order number	17.2112.201
Product description	Minivette® POCT Lithium heparin LH, nominal volume: 200 µl, preparation: Heparin, green, colour code ISO, 150 piece(s)/bag, 150 piece(s)/case

Product characteristics

Type of collection	capillary
Type of preparation	Heparin
Colour code	ISO
Application	capillary blood collection
Coding colour	green

Size

Sample volume	200 µl
---------------	--------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Minivette® POCT Lithium heparin LH, 200 µl,
plunger green, colour code ISO, 150 piece(s)/bag**

Page 2

Material & colours

Product material	Cellulose acetate (CA)
Colour of product	transparent

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Latex-free	yes
Batched	yes

Packaging

Minimum order qty.	1500
Type of smallest subpackaging	bag
Piece(s) / inner box	150
Piece (s) / outer case	1500
Piece(s) / pallet	108000
Depth of box	140 mm
Width of box	99 mm
Height of box	117 mm
Depth of case	513 mm
Width of case	293 mm
Height of case	140 mm
Case volume	0.021 cbm
Weight of product	0.0011 kg
Weight of case	2.37 kg
EAN of inner box	4038917339513
EAN case	4038917339506

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-04

This document was prepared by EDP support and is valid without signature.



SARSTEDT