



Product description

| | |
|---------------------|---|
| Order number | 17.2113.050 |
| Product description | Minivette® POCT EDTA K3E, nominal volume: 50 µl, preparation: K3 EDTA, red, colour code EU, 200 piece(s)/bag, 200 piece(s)/case |

Product characteristics

| | |
|---------------------|----------------------------|
| Type of collection | capillary |
| Type of preparation | K3 EDTA |
| Colour code | EU |
| Application | capillary blood collection |
| Coding colour | red |

Size

| | |
|---------------|-------|
| Sample volume | 50 µl |
|---------------|-------|

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Minivette® POCT EDTA K3E, 50 µl, plunger red,
colour code EU, 200 piece(s)/bag

Page 2

Material & colours

| | |
|-------------------|------------------------|
| Product material | Cellulose acetate (CA) |
| Colour of product | transparent |

Purity & certification

| | |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified | CE |
| Latex-free | yes |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 2000 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 200 |
| Piece (s) / outer case | 2000 |
| Piece(s) / pallet | 144000 |
| Depth of box | 140 mm |
| Width of box | 99 mm |
| Height of box | 117 mm |
| Depth of case | 513 mm |
| Width of case | 293 mm |
| Height of case | 140 mm |
| Case volume | 0.021 cbm |
| Weight of product | 0.001 kg |
| Weight of case | 2.68 kg |
| EAN of inner box | 4038917310406 |
| EAN case | 4038917310390 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-28

This document was prepared by EDP support and is valid without signature.



SARSTEDT