



Product Specification
End-to-end capillary, K2 EDTA, 20 µl, Ø: 1.28 mm, glass

Page 1

Product description

Order number	19.447.001
Product description	End-to-end capillary, K2 EDTA, capillary blood collection, preparation: K2 EDTA, nominal volume: 20 µl, Ø: 1.28 mm, glass, 100 piece(s)/container

Product characteristics

Type of collection	capillary
Type of preparation	K2 EDTA
Application	capillary blood collection
Volume of work	20 µl

Size

Sample volume	20 µl
Blood volume	20 µl
Diameter	1.28 mm
Length of product	29 mm

Material & colours

Product material	glass
Colour of product	transparent

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE supplier
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	container
Piece(s) / inner box	1000
Piece (s) / outer case	50000
Piece(s) / pallet	2400000
Depth of box	117 mm
Width of box	50 mm
Height of box	40 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
End-to-end capillary, K2 EDTA, 20 µl, Ø: 1.28 mm, glass

Depth of case	388 mm
Width of case	378 mm
Height of case	215 mm
Case volume	0.0315 cbm
Weight of product	0.0001 kg
Weight of case	8.65 kg
EAN of inner box	5706860174428
EAN case	5706860174435

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

