



Product description

Order number	20.1282
Product description	Microvette® 100 Lithium heparin LH, capillary blood collection, preparation: Lithium heparin LH, nominal volume: 100 µl, (LxØ) with cap: 47.6 x 10.8 mm, with print, cap orange, flat base, inner receptacle: conical, 100 piece(s)/case

Product characteristics

Type of collection	capillary
Type of preparation	Lithium heparin LH
Colour code	EU
Label/ Print	with print
Colour of print/label	orange
Cap	cap
Application	capillary blood collection
Base shape	flat base
Recommended centrifugation	2000 x g
Centrifugation time	5 Min
Shape of inner tube	conical

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
**Microvette® 100 Lithium heparin LH, 100 µl, cap
orange, flat base**

Page 2

Size

Sample volume	100 µl
Diameter	10.8 mm
Length including cap	47.6 mm
Length excluding cap	46.6 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	orange

Purity & certification

Product category	In vitro diagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	case
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	121000
Depth of box	139 mm
Width of box	108 mm
Height of box	66 mm
Depth of case	338 mm
Width of case	228 mm
Height of case	160 mm
Case volume	0.0123 cbm
Weight of product	0.0026 kg
Weight of case	3.11 kg
EAN of inner box	4038917086004
EAN case	4038917057134

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

