



## Product description

Order number	20.1292
Product description	Microvette® 200 Lithium heparin LH, capillary blood collection, preparation: Lithium heparin LH, nominal volume: 200 µl, (LxØ) with cap: 47.6 x 10.8 mm, with print, cap orange, flat base, inner receptacle: conical, 100 piece(s)/case

## Product characteristics

Type of collection	capillary
Type of preparation	Lithium heparin LH
Colour code	EU
Label/ Print	with print
Colour of print/label	orange
Cap	cap
Application	capillary blood collection
Base shape	flat base
Recommended centrifugation	2000 x g
Centrifugation time	5 Min
Shape of inner tube	conical

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Size

Sample volume	200 µl
Diameter	10.8 mm
Length including cap	47.6 mm
Length excluding cap	46.6 mm

## Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	orange

## Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

## Packaging

Minimum order qty.	1000
Type of smallest subpackaging	case
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	121000
Depth of box	139 mm
Width of box	108 mm
Height of box	66 mm
Depth of case	338 mm
Width of case	228 mm
Height of case	160 mm
Case volume	0.0123 cbm
Weight of product	0.0027 kg
Weight of case	3.23 kg
EAN of inner box	4038917086066
EAN case	4038917070386

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

