



**Product description**

|                     |  |
|---------------------|--|
| Order number        | 41.1438.901  |
| Product description | Micro sample tube, 500 µl Roche-Gen. 2/20 µl, capillary blood collection, preparation: Roche haemolysis gen. 2, incl. 20 µl end-to-end capillary, nominal volume: 20 µl, (LxØ) with cap: 42 x 10.8 mm, cap: white, push cap, skirted conical base, 2, 100 piece(s)/bag |

**Product characteristics**

|                     |  |
|---------------------|--|
| Type of collection  | capillary                                      |
| Type of preparation | Roche haemolysis gen. 2                        |
| Cap                 | flat cap                                       |
| Closure type        | push cap                                       |
| Application         | capillary blood collection                     |
| Example of use      | Blood sugar determination from capillary blood |
| Base shape          | skirted conical base                           |
| Volume of work      | 20 µl  |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Micro sample tube, 500 µl Roche-Gen. 2/20 µl, 20 µl, push cap**

**Page 2**

**Size**

|                      |         |
|----------------------|---------|
| Sample volume        | 20 µl   |
| Diameter             | 10.8 mm |
| Length including cap | 42 mm   |
| Length excluding cap | 39 mm   |

**Material & colours**

|                   |                                  |
|-------------------|----------------------------------|
| Product material  | Polypropylene (PP)               |
| Colour of product | transparent                      |
| Closure material  | Low Density Polyethylene (LD-PE) |
| Colour of cap     | white                            |

**Purity & certification**

|                  |                               |
|------------------|-------------------------------|
| Product category | In vitrodiagnostic device, CE |
| CE certified     | CE                            |
| Batched          | yes                           |

**Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 1000          |
| Type of smallest subpackaging | bag           |
| Piece(s) / inner box          | 100           |
| Piece (s) / outer case        | 1000          |
| Piece(s) / pallet             | 96000         |
| Depth of case                 | 348 mm        |
| Width of case                 | 278 mm        |
| Height of case                | 140 mm        |
| Case volume                   | 0.0135 cbm    |
| Weight of product             | 0.0021 kg     |
| Weight of case                | 2.42 kg       |
| EAN of inner box              | 4038917250269 |
| EAN case                      | 4038917250252 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

